

Compression Segments Head

General instructions

- Please bear in mind that the cutouts need to be kept as small as possible so as not to reduce the compression effect or therapy results. Too large cutouts could result in "window oedemas".
- In those styles with a closed forehead section, the entire back of the head is also closed.
In the case of the variant with an open forehead section, the back of the head is also open up to the height of the forehead.

Circumference measurements



Circumference "cA"

Lower neck circumference
parallel to the circumference
"cB"



Circumference "cB"

Neck circumference directly
below the chin and over the
first cervical vertebra at the
back of the head



Circumference "cC"

Over the tip of the chin and
over the first cervical vertebra
at the back of the head



Circumference "cD"

Directly below the lower lip
and over the first cervical ver-
tebra at the back of the head



Circumference "cE"

Directly above the eyebrows



Circumference "cE'"

Parallel to the circumference "cE"



Still Me Inc
Fitting Designs
Still Me Medical


225-273-1900 • Fax: 225-273-5555
info@stillmeinc.com

Length measurements



Length "LAB"

Width of the neck part: between "cB" and "cA", measuring at the front of the neck and following the body. Ensure that the distance is not so small that the fabric could roll up.



Length "LBC"

Directly below the chin to the tip of the chin, measuring at the front of the neck and following the body



Length "LCD"

From the tip of the chin to below the lower lip, measuring at the front and following the body



Length "LD'D1"

Headband from point "D1" over the largest area of the head to the point "D1" on the opposite side



Length "LEE"

Width of the headband: between "cE" and "cE1"



Length "LDE"

Position of the headband: between "cD" and "cE"

Additional measurements for face masks

Width measurements Please measure horizontally.

Length measurements Please measure vertically.



Width "K"

Mouth cutout, from one corner of the mouth to the other



Width "M"

Nose cutout, horizontally from the medial side of each nasal wing



Length "M2"

Root of the nose to the tip of the nose



Length "S"

Upper lip to the nose



Width "M1"

Knitted-in nose: from nasal wing to nasal wing over the tip of the nose measuring loosely



Width "N"

Eye cutout, from the corner of one eye to the corner of the other



Length "T"

Below the lower lip to the eye



Length "U"

Height of the eye cutout



Width "P"

Distance between the eyes



Width "Ear"

Ear cutout



Height "Ear"

Ear cutout