


## 225-273-1900 · Fax: 225-273-5555

info@stillmeinc.com

Patient: Last Name:	First Name
Fitter:	Clinic
Fitter Title:	_ (example: PT/OT/PTA)
Date:	



## Classic Arm Sleeve AG1 Ready-to-Wear





## Classic Arm Sleeve AG1 Ready-to-Wear

All measurements are in centimeters (cm).

Classic Arm Sleeve AG1 Size Chart					
	Least Wrist Circumference	Least Elbow Circumference Axilla Circumference			
	"C" Landmark	"E" Landmark	"G" Landmark		
Small	15.5–16.5	23-26	27.5-32		
Medium	16.5-18	26-29	32-36		
Large	18-19	29-33	36-41		
XLarge	19-20.5	33-37.5	41-45		

Classic Arm Sleeve AG1 Length Chart				
Wrist to Axilla (C to G) (measured with arm hanging down at side)				
Regular	43			
Long	45.5			









## Arion Easy-Slide Arm on 🖴

- The user-friendly application aid makes putting on compression arm sleeves quick and easy
- A straightforward donning method in combination with the application aid

Size	Circumference of widest part of the arm	BNR	UOM / Box	Order Qty.
Medium	14.5"-15.1" (37-38.5cm)	7966102	1	
Large	15.3"-16.1" (39-41cm)	7510001	1	