

Patient: Last Name: _____ First Name _____

Fitter: _____ Clinic _____

Fitter Title: _____ (example: PT/OT/PTA)

Date: _____

Compreflex® Reduce Calf

Previously **COMPREFLEX REDUCE BK**

sigvaris

Product Information

Product includes one Compreflex Reduce Calf and one pair of Cotton Liners.

☐ Right Calf Size: _____ Item #: _____
Color: ☐ Black ☐ Beige Length: _____ Quantity: _____

☐ Left Calf Size: _____ Item #: _____
Color: ☐ Black ☐ Beige Length: _____ Quantity: _____

Companion Products (sold separately)

☐ Compreflex Reduce Foot Color: ☐ Black ☐ Beige
Size: _____ Item #: _____ Quantity: _____
Size: _____ Item #: _____ Quantity: _____

Circumference

Left Right

C _____

A _____

I _____


J _____

G Length

Regular: 30–33.9cm
Tall: 34–42cm

G _____

K _____



K: Measure medial length from heel to 1st metatarsal head

Note: Compreboot® options available, each sold separately.

Measuring Instructions

Compreflex Reduce Calf:

Measure circumferences where indicated, beginning 5cm above the lateral malleolus (outer ankle bone), continuing to 25cm. Record in the corresponding lines A and C (Left/Right). Measure length from the lateral malleolus (outer ankle bone) to lateral aspect of the fibular head, and record in box G.

Compreflex Reduce Foot:

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

Compreflex Reduce Calf Size Chart (A–C)

	Range 1	Range 3
C	<65cm	<90cm
A	<55cm	<75cm

Black

Regular	2401-BKR	2403-BKR
Tall	2401-BKT	2403-BKT

Beige

Regular	2411-BKR	2413-BKR
Tall	2411-BKT	2413-BKT

Compreflex Reduce Foot Size Chart (I–K)

	Range 1	Range 3
I	30–44cm	39–54cm
J	24–33cm	30–43cm
K	16–25cm	20–27cm

Black

2801-BTL	2803-BTL
----------	----------

Beige

2811-BTL	2813-BTL
----------	----------