



Still Me Inc
Fitting Designs
Still Me Medical

225-273-1900 • Fax: 225-273-5555
info@stillmeinc.com

Patient: Last Name: _____ First Name _____

Fitter: _____ Clinic _____

Fitter Title: _____ (example: PT/OT/PTA)

Date: _____

sigvaris

Compreflex® Reduce Foot

Circumference


Left Right

I _____

J _____

Length

K _____



K: Measure medial length
from heel to 1st metatarsal head

Product Information

Product includes one Compreflex Reduce Foot.

☐ Right Foot

Size:

Item #:

Color: ☐ Black ☐ Beige

Quantity:

☐ Left Foot

Size:

Item #:

Color: ☐ Black ☐ Beige

Quantity:

Note: Compreflex Reduce Calf options available, each sold separately.

Compreflex Reduce Foot Size Chart (I-K)

	Range 1	Range 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
Black	2801-BTL	2803-BTL
Beige	2811-BTL	2813-BTL

Measuring Instructions

Compreflex Reduce Foot:

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.