

Patient: Last Name: _____ First Name _____

Fitter: _____ Clinic _____

Fitter Title: _____ (example: PT/OT/PTA)


Date: _____

sigvaris

Comprevest®

Step One


Measure the circumference in inches at the fullest part of bust.



Bust Size

Step Two

Measure the circumference in inches directly under the bust.



Band Size

Step Three

Determine the cup by subtracting the band from the bust.

Cup Size

2" or less = **A/B**
2" - 4" = **C/D**
4" - 6" = **DD/DDD**

Product Information

Product includes one Comprevest.

Comprevest Cup Size: Item #:
Color: Black White Band Size: Quantity:

Companion Products (sold separately)

Chip Pad Bra Half Item #:
Color: Black White Quantity:
 Chip Pad Bra Full Item #:
Color: Black White Quantity:
 Drain Tube Pouch Item #:
Color: Black White Quantity:

Caution: Elastic band contains natural rubber latex.

Comprevest Size Chart

	Small	Medium	Large	X Large	XX Large
Band Size	28-30in	32-34in	36-38in	40-42in	44-46in

Cup Size						
Black	A/B	2201-BRR-A	2202-BRR-A	2203-BRR-A	2204-BRR-A	2205-BRR-A
	C/D	2201-BRR-C	2202-BRR-C	2203-BRR-C	2204-BRR-C	2205-BRR-C
	DD/DDD	2201-BRR-D	2202-BRR-D	2203-BRR-D	2204-BRR-D	2205-BRR-D
White	A/B	2281-BRR-A	2282-BRR-A	2283-BRR-A	2284-BRR-A	2285-BRR-A
	C/D	2281-BRR-C	2282-BRR-C	2283-BRR-C	2284-BRR-C	2285-BRR-C
	DD/DDD	2281-BRR-D	2282-BRR-D	2283-BRR-D	2284-BRR-D	2285-BRR-D