



Still Me Inc
Fitting Designs
Still Me Medical

225-273-1900 • Fax: 225-273-5555
info@stillmeinc.com

Patient: Last Name: _____ First Name _____

Fitter: _____ Clinic _____

Fitter Title: _____ (example: PT/OT/PTA)

Date: _____


sigvaris

Neck and Mandible Pad

Circumference

Head _____

Neck _____



Product Information

Product includes one Neck and Mandible Pad.

☐ Neck & Mandible

Color: Black

Size:

Quantity:

Item #:

Neck & Mandible Pad Sizing Chart (Available in Black Only)

	Small	Medium	Large	X Large
Head	46–56cm	52–62cm	56–66cm	58–68cm
Neck	28–33cm	33–38cm	38–43cm	43–58cm
Black	2231-NM	2232-NM	2233-NM	2234-NM

Measuring Instructions

Head:

Measure the circumference by encircling the head; beginning at the top, moving below the chin; and finishing at the top of the head.

Neck:

Measure circumference at largest part of neck.