


225-273-1900 Fax: 225-273-5555 info@stillmeinc.com

SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)**
- 2 Calf Circumference at fullest part of the calf **(c)**
- 3 Calf Length from floor up to the bend of the knee (outside of calf) (A–D)

THIGH-HIGH STYLE

(continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference **(g)**
- 5 Thigh Length from floor up to the gluteal fold **(A–G)**
- 6 Hip Circumference at widest part of the hip **(f)**

Series	Features	Best Suited For		
SECURE BELOW KNEE	Maintains more graduated compression throughout the leg	Chronic venous insufficiency (CVI)		
78% Nylon	compared to traditional stockings	Lymphedema or Lipedema		
22% Spandex	High modulus in-lay yarns provide containment and work to maintain	Recurrent venous ulcers		
(contents may very by style)	limb profile	Deep vein thrombosis (DVT)		
Not made with	Secure dotted grip-top band Double-covered yarns ensure no	Superficial phlebitis		
natural Latex rubber	bare spandex next to the skin	Post-thrombotic syndrome		
	Unique sizing system fits a wide	Severe varicosities		
	range of regular and full sizes— up to 40 inch thighs	Orthostatic hypotension and post-phlebotic syndrome		

Donning Accessories				
	SIGVARIS ARION Sim-Slide (Open Toe Donning Device)		SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device)	
	Latex-Free Donning Gloves		SIMON (Donning Device)	
	Rubber Donning Gloves		SIGVARIS Doff N' Donner (Donning Device)	
	Cotton Underliner		S.O.S. (Slip On SIGVARIS Donning Device	

Products	Compression	Style		Colors
550 SECURE BELOW KNEE*	20-30 mmHg 30-40 mmHg 40-50 mmHg	Closed Toe: ☐ Calf w/grip-top ☐ Thigh w/grip-top	☐ Beige (77)	□ Black (99)

 $^{^{\}star}$ 40–50mmHg made on demand. Please allow 7–10 days for delivery.